Hindawi Publishing Corporation Journal of Applied Mathematics Volume 2012, Article ID 591252, 12 pages doi:10.1155/2012/591252

# Research Article

# Monitoring Personalized Trait Using Oscillometric Arterial Blood Pressure Measurements

# Young-Suk Shin

Department of Information and Communication Engineering, Chosun University, 375 Seosuk-dong, Dong-gu, Gwangju 501-759, Republic of Korea

Correspondence should be addressed to Young-Suk Shin, ysshin@chosun.ac.kr

Received 1 October 2011; Accepted 29 November 2011

Academic Editor: Pedro Serranho

Copyright © 2012 Young-Suk Shin. This is an open access article distributed under the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

The blood pressure patterns obtained from a linearly or stepwise deflating cuff exhibit personalized traits, such as fairly uniform peak patterns and regular beat geometry; it can support the diagnosis and monitoring of hypertensive patients with reduced sensitivity to fluctuations in Blood Pressure (BP) over time. Monitoring of personalized trait in Oscillometric Arterial Blood Pressure Measurements (OABPM) uses the Linear Discriminant Analysis (LDA) algorithm. The representation of personalized traits with features from the oscillometric waveforms using LDA algorithm includes four phases. Data collection consists of blood pressure data using auscultatory measurements and pressure oscillations data obtained from the oscillometric method. Preprocessing involves the normalization of various sized oscillometric waveforms to a uniform size. Feature extraction involves the use of features from oscillometric amplitudes, and trait identification involves the use of the LDA algorithm. In this paper, it presents a novel OABPM-based blood pressure monitoring system that can monitor personalized blood pressure pattern. Our approach can reduce sensitivity to fluctuations in blood pressure with the features extracted from the whole area in oscillometric arterial blood pressure measurement. Therefore this technique offers reliable blood pressure patterns. This study provides a cornerstone for the diagnosis and management of hypertension in the foreseeable future.

## 1. Introduction

Blood pressure (BP) is a vital sign, which along with body temperature, heart rate, and respiratory rate provides various physiological statistics about the body. Small changes in the BP over a period of time can provide clues about cardiovascular and respiratory abnormalities in a patient. Oscillometry is one of the widely used methods to determine the blood pressure [1–7]. The oscillometric method of measuring blood pressure uses the amplitude of cuff pressure oscillations from a linearly or stepwise deflating cuff and is given as two values, the systolic and diastolic pressures. The cuff pressure oscillations consist of waveforms. The systolic

pressure is the pressure associated with contraction of the heart, and indicates the maximum amount of work per stroke needed for the heart to pump blood through the arteries [8]. In contrast, diastolic pressure is the pressure in the large arteries during relaxation of the heart left ventricle [9]. The diastolic pressure indicates the amount of pressure that the heart must overcome in order to generate the next beat [8].

There have been ongoing studies to develop reliable measurements of blood pressure [3–7]. These researches have focused on improving the accuracy of blood pressure measurements. However, a large number of cardiovascular diseases such as arrhythmia can make it difficult to obtain accurate blood pressure measurements [3]. To determine the true BP level, many BP measurements need to be taken over a long period of time and problems affected by the white-coat effect have to be solved. The white-coat effect is usually defined as the difference between the BP measured at home and at the office. White-coat effect can be influenced by anxiety, a hyperactive alerting response, or a conditioned response. The white-coat effect typically causes the office BP to be higher than the home BP and is present in a high percentage of hypertensive patients [10]. If there are personalized traits in blood pressure measurements, problems such as noises caused by cardiovascular diseases like arrhythmia or problems of the white-coat effect may be overcome. Therefore, this study proposes the oscillometric measurement-based automatic blood pressure pattern identification system to explore personalized traits prior to obtaining reliable blood pressure measurements. The proposed approach demonstrates the feasibility of personalized trait identification with 85 people.

This paper aims to explore blood pressure pattern identification to find personalized traits in oscillometric arterial blood pressure measurements using the linear discriminant Analysis (LDA) algorithm. Section 2 introduces a review of related work. Section 3 develops a representation of personalized traits with features from the oscillometric waveforms. It consists of four steps. The first step introduces the database used for this research. The second step presents a preprocessing technique for obtaining uniform sized oscillation waves, and the third step develops a personalized traits representation via oscillations of amplitude features from uniform sized oscillation waves. The fourth step describes data reduction and feature extraction using LDA in the appearance-based approach. Section 4 presents the performance of the blood pressure patterns identification model via the LDA algorithm. Finally, this study discusses the advantages and applications of personalized trait monitoring.

## 2. Related Work

Blood pressure best predicts cardiovascular risk. Therefore, a variety of studies have been proposed to improve the accuracy of blood pressure measurements [3–7, 11–14]. Many studies use the oscillometric method to measure the blood pressure [1–7, 11, 12]. The oscillometric method is used to find the peak values of the oscillation waveform, which are determined as the oscillation amplitudes obtained from the pressure of the linearly deflating cuff. This method has virtually no complications and needs less expertise; it is less unpleasant and painful for the patient. In [1–7, 11, 12], blood pressure measurements based on the oscillometric method typically only use single-point estimates for both systolic blood pressure and diastolic blood pressure. Recently, BP measurements in [13, 14] were introduced: the confidence interval estimate of the systolic blood pressure and diastolic blood pressure. In [13], the confidence interval estimate performed well only when sample size is large. The confidence interval estimate used in [14] requires independent and identically distribution of data. But these methods also have to measure single-point estimates for systolic and diastolic blood pressure and can reflect on sensitivity to fluctuations in BP measurements.

In this paper, we have attempted to extract personalized blood pressure patterns of oscillation amplitudes rather than measure single-point estimates for systolic and diastolic blood pressure. During feature extraction, we focus on the more uniform features of the oscillation amplitudes in each person.

# 3. Methodology

This section describes a new blood pressure patterns identification technique to find personalized traits in oscillometric arterial blood pressure measurements using the LDA algorithm. This work consists of four steps. First, data collection is described. Second, oscillometric waveforms of various sizes are normalized to a uniform size. Third, features based on the oscillation amplitudes are developed. Finally, the LDA algorithm is applied to identify blood pressure patterns.

#### 3.1. Data Collection

Experimental data has been provided by the blood pressure research team of the University of Ottawa in Canada. The database consists of blood pressure data using auscultatory measurements and pressure oscillations data obtained from the oscillometric method. The blood pressure data measured using the auscultatory method was obtained by two trained nurses. The oscillometric method is similar to the auscultatory technique, but it uses a pressure sensor instead of a stethoscope to record the pressure oscillations within the cuff. This method requires an external inflatable cuff, which can be placed around the left wrist at heart level. The cuff is inflated starting from below the diastolic pressure until the cuff pressure exceeds the systolic pressure. The cuff pressure is first increased until it exceeds the systolic pressure and then deflated until it reaches certain pressures at fixed or variable intervals [7].

The database consists of a total of  $425~(85\times5)$  records with five recordings per subject from 85 male and female subjects. Subjects met various blood pressure criteria: 10% of participants had BP below 100 mmHg systolic, 10% had BP above 140 mmHg systolic, 10% had BP below 60 mmHg diastolic, 10% had BP above 100 mmHg diastolic, and the remainder had BP distributed between these outer limits. The subjects' ages ranged from 10 to 80 years. Subjects were allowed to relax in a waiting room area for 15 minutes and the measurement room was organized to be conducive to accurate blood pressure measurements. The subjects were told not to talk or move during the readings. Five records per subject were acquired, and measurements were repeated for one minute with a one-minute rest period. Figure 1 shows one example of an oscillation pattern extracted from the cuff pressure acquired from the oscillometric method.

#### 3.2. Preprocessing

The number of oscillation waveforms extracted from the cuff pressure varies according to physiology, geometry of the heart, hypertension, gender, and age (see Table 1). Table 1 shows a partial example of varying number of oscillation waveforms extracted from the cuff pressure. The systolic and diastolic pressures are the average values acquired by two nurses with auscultatory measurements. We can find a similar number of oscillation waveforms in 5 measurements of the same subject. That is, the same person can have similar number of oscillation waveforms. We attempt to use normalization to reduce variations of corresponding




Figure 1: Oscillation pattern extracted from the cuff pressure with the oscillometric method.

oscillation waveforms for different oscillation waveforms of the same person. A blood pressure pattern means a varying number of oscillation waveforms in one record for one-minute. Training set is defined as X. Given the training set  $X = \{W_i\}_{i=1}^N$ , containing N blood pressure patterns where each blood pressure pattern  $W_i = \{W_{ij}\}_{j=1}^{N_i}$  consists of a number of oscillation waveforms  $W_{ij}$ , the normalization is applied as follows:

$$\varphi = \operatorname{sqrt}\left(\sum_{j=1}^{N_i} (W_{ij})^2\right),$$

$$W_i^* = \frac{W_i}{\varphi}.$$
(3.1)

## 3.3. Feature Extraction

The proposed feature extraction technique extracts features of mean amplitude (MA), maximum positive amplitude (MPA), and maximum negative amplitude (MNA) based on database with the number of oscillation waveforms.

To implement the proposed approach, we segment a normalized oscillation pattern into 29-sample windows (at least including a single beat in the minimum oscillation waveforms) to obtain the feature windows. That is, a blood pressure pattern is divided into 29 sections and each divided section has to include at least a single heartbeat. In this study, 29 sections are defined for including at least a single heart beat on the training set, *X*. Blood pressure patterns larger than the minimum number of oscillation waveforms in a blood pressure pattern represent multiple heart beats within a given window. One-feature window means one section in 29 sections. Figure 2 shows four heart beats detected within a given

Table 1: Oscillation waveforms of various sizes extracted from the cuff pressure.

S1         50         Male         2         4679         128         98           S1         5         4679         128         98         Prehypertensic           S2         4         4076         126         98         Prehypertensic           S2         1         4         4075         126         98         Prehypertensic           S2         2         4684         125         99         Prehypertensic           S2         2         3499         114         68         Normal           S3         54         4         3523         142         88         Prehypertensic           S3         54         4         362         145         88         Prehypertensic           S4         34         345         145         88         Prehypertensic           S4         34         345         145         88         Prehypertensic           S4         34         340         87         Prehypertensic           S4         34         848         131         82           S5         36         4         824         131         134         134           S6 <t< th=""><th>Subjects</th><th>age</th><th>Gender</th><th>Reading number</th><th>Number of oscillation waveforms acquired from cuff pressure for oscillometric method</th><th>Blood pressure by auscultatory method Systolic Diastolic</th><th>uscultatory method Diastolic</th><th>Note</th></t<>	Subjects	age	Gender	Reading number	Number of oscillation waveforms acquired from cuff pressure for oscillometric method	Blood pressure by auscultatory method Systolic Diastolic	uscultatory method Diastolic	Note
50         Male         2         4679         128         98           4         4         4084         126         98           5         4684         126         98           6         4         4085         126         98           1         3457         113         67         99           2         3457         113         68         68           1         4         3457         104         68           1         5         3551         104         66           1         2         3845         142         88           1         34         382         142         88           1         3879         142         88           1         8401         131         86           2         8500         131         86           3         850         130         77           4         8808         131         78           5         4         808         131         78           6         4         4250         109         69           7         4         4220         109				1	4591	127	86	
50         Male         3         4684         126         98           4         4         4065         126         98           5         4         4698         126         98           1         1         3457         113         67           2         3499         104         68           4         3551         104         66           5         3499         104         68           6         3         3551         104         61           8         2         3552         104         61           9         3         3553         142         88           1         2         3842         142         88           1         34         3892         143         88           1         4         3955         140         89           1         8500         131         86           1         8500         131         86           1         804         130         17           1         804         130         104         77           1         4         4220         104				2	4679	128	86	
22         Female         3         3457         1126         98           22         1         3457         113         67           24         3459         104         68           25         3450         104         66           25         3551         104         66           25         3551         104         66           26         3551         104         67           27         3721         145         88           28         3845         142         88           34         Male         3         395         143         86           34         Male         3         874         133         86           34         Male         3         874         133         78           34         Male         3         874         133         79           4         4         8808         131         78           5         4         4069         106         71           6         4         4109         104         72           7         4         4106         109         69	S1	20	Male	3	4684	126	86	Prehypertension
22         Female         3         4598         125         99           22         3497         113         67         68           3         3530         1104         68         68           4         3551         104         61         66         68           54         1         3721         104         61         88         65         65         65         65         66         67				4	4705	126	86	4
22         Female         2         3457         113         67           22         Female         3         3330         104         68           4         3551         104         66         66           54         3552         108         65         88           54         5         3845         145         88           54         845         142         88         88           54         845         140         88         88           34         Male         3         8746         130         86           34         Male         3         8746         130         77           34         Male         3         8746         130         77           34         Male         3         8746         130         77           35         Male         3         8746         130         77           36         Male         3         4069         106         71           36         Male         4         4208         108         70           43         Female         3         4041         109         69				Ŋ	4698	125	66	
22         Female         2         3499         104         68           4         3550         112         66           5         3551         104         61           5         3551         104         66           1         3521         104         65           2         3845         142         88           4         3892         143         84           5         3979         143         80           34         Male         2         8401         131         86           34         4         8808         131         86           34         4         8808         131         77           36         44         8808         131         78           36         4         8948         130         71           36         4         8948         130         71           36         4         4069         104         72           36         44520         109         67           4         4220         109         69           4         4220         112         69				₽	3457	113	29	
22         Female         3         3530         112         66           4         3551         104         61         61           5         3551         104         61         61           6         3         3523         108         65           54         3845         142         88         83           1         4         3892         144         88           1         4         3892         144         88           1         3797         131         82           1         8401         131         82           8500         131         82           8500         131         82           846         131         73           850         131         73           850         131         73           848         130         82           848         130         67           848         4459         104         71           85         4409         108         70           85         4410         109         69           85         4410         112         69     <					3499	104	89	
54       3551       104       61         54       3523       108       65         1       3721       145       88         2       3845       143       88         4       3935       143       84         5       3879       140       80         34       Male       3       870       131       86         4       8808       131       82         5       8948       130       79         6       8948       130       77         7       1       3857       104       72         8       4       4220       104       72         4       4220       108       70         4       4220       108       70         4       4220       108       70         4       4220       108       69         4       430       4413       109       70         4       430       440       109       70         8       69       112       69         9       44139       114       71         11       114       71	S2	22	Female		3530	112	99	Normal
54         Female         1         3721         145         88           54         3845         142         88         83           54         3892         142         88           5         3892         140         84           1         4         3892         140         84           2         8401         135         78           34         8401         131         86           5         8500         131         86           1         8746         130         79           1         8746         130         79           1         8746         130         77           808         131         86           808         131         86           808         130         77           808         130         70           808         109         77           809         104         72           80         4         4220         109         70           81         106         70         109         70           82         4220         108         70         70 <td></td> <td></td> <td></td> <td>4</td> <td>3551</td> <td>104</td> <td>61</td> <td></td>				4	3551	104	61	
54     Female     1     3721     145     88       3845     142     83       4     3845     142     83       5     3892     143     84       1     8401     135     78       2     8500     131     86       34     Male     3     8746     130     79       5     8808     130     87       8808     130     87       1     3957     109     67       1     3957     104     72       4     4220     103     71       4     4220     103     71       4     4220     103     69       43     Female     3     404     112     69       4     4139     410     114     71       69     114     71     69				Ŋ	3523	108	65	
54         Female         2         3845         142         83           4         3892         143         84         84           5         3935         140         84         84         84         84         84         84         84         84         84         86					3721	145	88	
54         Female         3         3892         143         84           4         3935         140         80         80           34         Male         2         8401         131         82           34         Male         3         8746         131         86           34         808         131         78         82           34         1         808         131         78           36         1         3957         109         67           4         4069         106         71           4         420         104         72           4         420         103         70           4         420         103         70           4         4404         109         70           4         4404         109         70           4         4139         112         69           4         4139         114         71				2	3845	142	83	C40.00 1
34       3935       140       80         35       3979       135       78         34       Male       3       8746       131       86         4       8808       131       78         8948       131       78         36       Male       3       4159       67         4       4220       104       72         4       4220       103       71         43       Female       3       4106       103       70         43       Female       3       4041       109       69         43       Female       3       4041       109       69         43       Female       3       4041       109       69         43       Female       3       4106       109       69         43       Female       3       4041       109       70         43       Female       3       4041       109       69         44       4139       69       71       69         44       4139       69       71       69         44       4139       69       71       71	S3	54	Female		3892	143	84	Stage 1
34       Male       3       8401       131       82         34       Male       3       8746       130       79         4       8808       131       82         5       8948       130       79         1       3957       109       67         2       4069       104       72         4       4220       103       71         5       4218       108       70         43       Female       3       4041       109       69         43       Female       3       4041       109       69         43       Female       3       4041       109       69         43       5       4041       112       69         44       4139       112       69         5       4041       110       69				4	3935	140	80	nypertension
34     Male     3     8401     131     82       5     8500     131     86       5     808     131     78       5     8948     130     79       1     3957     109     67       2     4069     106     71       4     4220     103     71       5     4218     108     70       43     Female     3     4041     109     69       43     Female     3     4041     112     69       4     4139     114     71				Ŋ	3979	135	78	
34       Male       3       8500       131       86         4       8808       131       78         5       8948       130       78         1       3957       109       67         2       4069       106       71         36       Male       3       4159       104       72         4       4220       103       71         5       4218       108       70         43       Female       3       4041       109       69         43       Female       3       4041       112       69         4       4139       4139       114       71				1	8401	131	82	
34     Male     3     8746     130     79       4     8808     131     78       5     8948     130     78       1     3957     109     67       2     4069     106     71       4     4220     103     71       5     4218     108     70       43     Female     3     4041     109     69       43     Female     3     4139     112     69       43     44     4139     114     71				2	8500	131	98	
4     8808     131     78       5     8948     130     82       1     3957     109     67       2     4069     106     71       4     4220     103     72       5     4218     108     70       43     Female     3     4041     109     69       43     Female     3     4041     109     70       5     4139     64     70       44     4139     112     69       5     4100     114     71	S4	34	Male	3	8746	130	26	Prehypertensive
36     Male     3     4069     67       36     Male     3     4159     104     72       43     4220     103     71       5     4218     108     70       43     Female     3     4041     109     69       43     Female     3     4041     109     70       4     4139     64     70       4     4139     112     69       5     4100     114     71				4	8808	131	78	
36     Male     3     4069     67       36     Male     3     4159     104     72       4     4220     103     71       5     4218     108     70       1     3991     109     69       43     Female     3     4041     109     70       4     4139     641     71       5     4100     112     69       7     4139     114     71				rC	8948	130	82	
36     Male     3     4069     106     71       4     4     4220     104     72       5     4218     108     70       43     Female     3     4041     109     69       4     4139     112     69       4     4139     114     71					3957	109	29	
36     Male     3     4159     104     72       4     4220     103     71       5     4218     108     70       1     3991     109     69       2     4106     112     69       4     4139     112     69       5     4100     114     71				2	4069	106	71	
4     4220     103     71       5     4218     108     70       1     3991     109     69       2     4106     112     69       4     4139     112     69       5     4100     114     71	S <sub>5</sub>	36	Male	3	4159	104	72	Normal
4218     108     70       1     3991     109     69       2     4106     112     69       4     4041     109     70       4     4139     112     69       5     4100     114     71				4	4220	103	71	
1     3991     109     69       2     4106     112     69       4     4041     109     70       4     4139     112     69       5     4100     114     71				Ŋ	4218	108	70	
43 Female 3 4106 112 69 4 4041 109 70 4 4139 112 69 5 4100 114 71				П	3991	109	69	
43 Female 3 4041 109 70 4 4139 112 69 5 4100 114 71				2	4106	112	69	
4139 112 4100 114	9S	43	Female		4041	109	70	Normal
4100 114				4	4139	112	69	
				5	4100	114	71	





Figure 2: Four heartbeats detected in the segmented one-feature window from one subject.

window. We extract three features from the oscillation amplitudes in the segmented feature window: mean amplitude, maximum positive amplitude, and maximum negative amplitude. The circle marks of Figure 2 display the maximum positive and negative amplitudes extracted in the segmented feature window. We define the following for three features; MA means the averaged oscillation amplitudes in the segmented one-feature window; MPA means the amplitude of the maximum upper pulse from the oscillations in the segmented one-feature window; MNA describes the amplitude of the maximum lower pulse from the oscillations in the segmented one-feature window.


Figure 3 shows the feature extraction results of one subject with the mean amplitude and maximum positive and negative amplitudes in each feature window. The mean amplitude can reduce noise signals within the feature window, and the maximum positive and negative amplitudes exhibit personalized traits in the period of high or low cuff pressure. Figures 4 and 5 show the feature extraction results of six subjects in Table 1 with maximum positive and negative amplitudes in 29 feature windows, respectively. Figures 4 and 5 show the averaged results of five readings obtained from the oscillometric blood pressure measurements of each subject. Subjects S2, S5, and S6 of Figures 4 and 5 are normal (BP: <120/80 mmHg), whereas blood pressure subjects S1 and S4 are prehypertensive (BP: 120/80 to 139/89 mmHg) and S3 is stage 1 hypertensive (BP: 140/90 to 160/100 mmHg) blood pressure subject. In Figures 4 and 5, stage 1 hypertensive or prehypertensive subjects display a steep-slope pattern in front of the feature windows compared to normal subjects. Especially, older subjects show higher amplitudes based on the MPA features. In the MNA features, stage 1 hypertensive or prehypertensive or prehypertensive compared to normal subjects.

## 3.4. Identification


Linear discriminant analysis is used for data reduction and feature extraction in the appearance-based approach. LDA searches for feature vectors in the fundamental space that best discriminates among classes [15]. LDA describes a linear combination of feature vectors that produces the largest mean differences between the target classes. Features of the heartbeat applied for human identification from electrocardiogram (ECG) [16–18] are similar to features extracted from oscillometric arterial blood pressure measurements. Features extracted from oscillometric measurements are classified with an appearance-based approach based on LDA. Appearance-based approach is usually taken by different two-dimensional views



**Figure 3:** Features extracted with mean amplitude and maximum positive and negative amplitudes in the segmented feature windows from one subject.



**Figure 4:** Features extracted with maximum positive amplitude averaged from five readings in the feature windows for normal and hypertension (prehypertension) subjects with respect to age.



**Figure 5:** Features extracted with maximum negative amplitude averaged from five readings in the feature windows for normal and hypertension (prehypertension) subjects with respect to age.

of the object of interest. These methods based on the applied features can be subdivided into two approaches: local and global approaches. This study applies global appearance-based method. The main idea is to project the original input data onto a suitable lower-dimensional subspace that represents the data best for a specific work. Selecting optimization criteria for the projected data is the goal to best identify personalized trait.

Given a training set  $X = \{W_i^*\}_{i=1}^C$ , containing C classes with each class  $W_i^* = \{w_{ij}^*\}_{j=1}^{C_i}$  consisting of a number of features,  $w_{ij}^*$ , there are a total of  $N = \sum_{i=1}^{C} C_i$  oscillation patterns. We define two measures for all samples of all classes.  $S_{WT}$  is defined as within-class scatter matrices of the training feature set.  $S_{BT}$  is defined as between-class scatter matrices of the training feature set.  $S_{WT}$  and  $S_{BT}$  are given as

$$S_{WT} = \frac{1}{N} \sum_{i=1}^{c} \sum_{j=1}^{C_i} (w_{ij}^* - \mu_i) (w_{ij}^* - \mu_i)^T,$$

$$S_{BT} = \frac{1}{N} \sum_{i=1}^{c} (\mu_i - \mu) (\mu_i - \mu)^T.$$
(3.2)

In (3.2),  $w_{ij}^*$  denotes the jth sample of class i, c is the number of classes,  $\mu_i$  is the mean of class i, and  $c_i$  denotes the number of samples in class i and  $\mu$  is the mean of all classes. The LDA approach [19] finds a set of basis vectors described by  $\varphi$  that maximizes the ratio between  $S_{BT}$  and  $S_{WT}$ :

$$\varphi = \arg\max \frac{|\varphi^T S_{\text{BT}} \varphi|}{|\varphi^T S_{\text{WT}} \varphi|}.$$
(3.3)

One method is to assume that  $S_{\rm WT}$  is nonsingular and the basis vectors  $\varphi$  correspond to the first N eigenvectors with the largest eigenvalues of  $S_{\rm WT}^{-1}S_{\rm BT}$ . LDA-based feature representation,  $y = \varphi^T w^*$ , is produced by projecting the normalized input features ( $w^*$ ) from the oscillation amplitudes onto the subspace spanned by the N eigenvectors.


# 4. Experimental Results

To evaluate the performance of our approach, we conducted our experiments with the pressure oscillations data (measured using oscillometric method) provided by the blood pressure research team of the University of Ottawa in Canada. For the experiment, we used 425 records with five readings per subject obtained from 85 subjects; the training set consisted of 255 records with three readings per subject obtained from 85 subjects; the testing set consisted of the remaining 170 data readings excluded from the training set, two readings per subject obtained from 85 subjects. The blood pressure data measured using the auscultatory method provided indirect information for analyzing the subjects recognized by the oscillometric method in our experiment.

The experimental results were evaluated with the performance of LDA by using the nearest neighbor algorithm. The Euclidean distance was used for the similarity measure. To find the optimal LDA-based features, our implementation used the five sets of features from Figure 3 to test their discrimination power. One set included all of the features, whereas

Subset	Feature
I	Mean amplitude (MA)
II	Maximum positive amplitude (MPA)
III	Maximum negative amplitude (MNA)
IV	MPA + MNA
V	MA + MPA + MNA

**Table 2:** Subsets of features extracted using the oscillometric method.



**Figure 6:** Features extracted with mean amplitude and maximum positive and negative amplitudes in the feature windows from normal and hypertension (prehypertension) subjects.

the other four sets included a subset of these, as shown in Table 2. Subset IV concatenates the features of the maximum positive and negative amplitudes, and subset V concatenates the features of the mean amplitudes and maximum positive and negative amplitudes into one vector.

Figure 6 shows the feature extraction results obtained from six subjects in Table 1 and the mean amplitudes and maximum positive and negative amplitudes are concatenated into one vector. S2, S5, and S6 are normal blood pressure subjects, and S1, S3, and S4 are prehypertension or hypertensive blood pressure subjects. This shows the averaged results of the five readings for the oscillometric blood pressure measurements obtained from each subject. The feature windows describe feature windows 1 to 29 extracted from the maximum positive amplitudes, feature sections 30 to 57 extracted from the maximum negative amplitudes, and feature sections 58 to 87 extracted from the mean amplitude. The stage 1 hypertensive subject (S3) displays steeper maximum positive amplitude than that of the normal subjects. Prehypertensive or stage 1 hypertensive subjects generally display lower maximum negative amplitude than that of normal subjects. This shows that the averaged features of the five readings taken from each subject are plotted in a personalized uniform pattern. The results of the final LDA-based experiments are listed in Table 3. We can see that using all of the features provides the best blood pressure pattern identification rate, and subset IV shows good performance, while subset I shows the worst performance. LDA does not go beyond 85 for the dimensionality of the LDA space. Since we use 85 classes, this gives us an upper bound of 85-dimensional LDA space.

Subset	Recognition rate (%)
I	34.30
II	67.44
III	72.09
IV	93.02
V	94.70

**Table 3:** Experimental results of LDA.




Figure 7: LDA recognition performance according to dimensionality via LDA algorithm with nearest neighbor classifier.

Figure 7 shows the recognition results based on the dimensionality that yields the best identification rate. We achieved the best blood pressure pattern identification rate of 94.7% for the first 18 eigenvectors. Thus, the first 18 eigenvectors are estimated to the optimal decision boundary to best identify personalized trait using LDA in this study.

#### 5. Discussion

This study aimed to explore a new blood pressure patterns identification model for personalized traits monitoring of oscillometric arterial blood pressure measurements using the linear discriminant analysis algorithm. A blood pressure patterns identification model was used for the oscillometric arterial blood pressure measurements, which successfully discriminated personalized traits for the LDA algorithm. Our best recognition result showed a recognition rate of 94.7% for the first 18 eigenvectors. This means that the optimal LDA-based 18 eigenvectors in oscillometric arterial blood pressure measurements can effectively represent personalized traits.

The personalized traits of the oscillometric arterial blood pressure measurements can be represented for the features extracted from the whole domain of one oscillation pattern. Especially, the integration of the three feature streams extracted from each segmented feature window for the whole domain of one oscillation pattern enhances the recognition performance. In our experiment, the integration of the feature streams extracted with the maximum positive and negative amplitudes largely improved the recognition rate. In the three feature streams, while the maximum positive and negative amplitude feature streams showed strong effects on the recognition performance, the mean amplitude showed a weak effect. We propose that the maximum positive and negative amplitude features can effectively represent personalized traits of oscillometric arterial blood pressure measurement.

Features extracted from each segmented feature window in the oscillometric method may support the monitoring and diagnosis of hypertensive patients because stage 1 hypertensive or prehypertensive subjects display a steep-slope pattern in front of the feature windows compared to normal subjects.

Our approach offers a simple and inexpensive means of monitoring personalized trait with blood pressure patterns in oscillometric arterial blood pressure measurement. Based on these results, this study has established a new blood pressure monitoring system for health care monitoring in oscillometric arterial blood pressure measurements. Our research has the potentiality for the diagnosis and management of hypertension and provides a foundation of a new biometric modality using blood pressure patterns.

# Acknowledgments

This study was supported by research fund from Chosun University, 2011. The author wishes to acknowledge the data support of the blood pressure research team of the University of Ottawa in Canada.


#### References


- [1] L. A. Geddes, M. Voelz, and C. Combs, "Characterization of the oscillometric method for measuring indirect blood pressure," *Annals of Biomedical Engineering*, vol. 10, no. 6, pp. 271–280, 1982.
- [2] J. C. T. B. Moraes, M. Cerulli, and P. S. Ng, "Development of a new oscillometric blood pressure measurement system," in *Proceedings of the 26th IEEE Computers in Cardiology Conference*, pp. 467–470, September 1999.
- [3] C. T. Lin, S. H. Liu, J. J. Wang, and Z. C. Wen, "Reduction of interference in osszcillometric arterial blood pressure measurement using fuzzy logic," *IEEE Transactions on Biomedical Engineering*, vol. 50, no. 4, pp. 432–441, 2003.
- [4] T. J. Dorsett, "Application of a prediction and smoothing algorithm to non-invasive blood pressure measurement," in *Proceedings of the 13th Annual International Conference of the IEEE Engineering in Medicine and Biology Society*, pp. 468–469, November 1991.
- [5] K. G. Ng and C. F. Small, "Survey of automated noninvasive blood pressure monitors," *Journal of Clinical Engineering*, vol. 19, no. 6, pp. 452–475, 1994.
- [6] T. J. Brinton, B. Cotter, M. T. Kailasam et al., "Development and validation of a noninvasive method to determine arterial pressure and vascular compliance," *American Journal of Cardiology*, vol. 80, no. 3, pp. 323–330, 1997.
- [7] S. Colak and C. Isik, "Blood pressure estimation using neural networks," in *Proceedings of the IEEE International Conference on Computational Intelligence for Measurement Systems and Applications*, pp. 14–16, Boston, Mass, USA, 2007.
- [8] E. P. Dutton, "Heartbook, a guide to prevention and treatment of cardiovascular disease. The American heart association," vol. 116, pp. 132–136, 1980.
- [9] N. Kaplan, Clinical Hypertension, Williams and Wilkins, Baltimore, Md, USA, 3nd edition, 2005.
- [10] T. G. Pickering, N. H. Miller, G. Ogedegbe, L. R. Krakoff, N. T. Artinian, and D. Goff, "Call to action on use and reimbursement for home blood pressure monitoring: a joint scientific statement from the American heart association, american society of hypertension, and preventive cardiovascular nurses association," *Journal of Cardiovascular Nursing*, vol. 23, no. 4, pp. 299–323, 2008.
- [11] M. Forouzanfar, H. R. Dajani, V. Z. Groza, M. Bolic, and S. Rajan, "Oscillometric blood pressure estimation using principal component analysis and neural networks," in 2009 IEEE Toronto International Conference and Technology for Humanity (TIC-STH '09), pp. 981–986, September 2009.
- [12] S. Chen, V. Z. Groza, M. Bolic, and H. R. Dajani, "Assessment of algorithms for oscillometric blood pressure measurement," in *Proceedings of the IEEE Intrumentation and Measurement Technology Con*ference, pp. 1763–1767, May 2009.
- [13] L. R. Krakoff, "Confidence limits for interpretation of home blood pressure recordings," Blood Pressure Monitoring, vol. 14, no. 4, pp. 172–177, 2009.

- [14] S. Lee, M. Bolic, V. Z. Groza, H. R. Dajani, and S. Rajan, "Confidence interval estimation for oscillometric blood pressure measurements using bootstrap approaches," *IEEE Transactions on Instrumentation and Measurement*, vol. 60, no. 10, pp. 3405–3415, 2011.
- [15] K. Fukunaga, Introduction to Statistical Pattern Recognition, Academic Press, 2nd edition, 1990.
- [16] Y. Wang, F. Agrafioti, D. Hatzinakos, and K. N. Plataniotis, "Analysis of human electrocardiogram for biometric recognition," EURASIP Journal on Advances in Signal Processing, vol. 2008, Article ID 148658, 11 pages, 2008.
- [17] S. A. Israel, J. M. Irvine, A. Cheng, M. D. Wiederhold, and B. K. Wiederhold, "ECG to identify individuals," *Pattern Recognition*, vol. 38, no. 1, pp. 133–142, 2005.
- [18] C. C. Chiu, C. M. Chuang, and C. Y. Hsu, "A novel personal identity verification approach using a discrete wavelet transform of the ECG signal," in *Proceedings of the International Conference on Multi*media and Ubiquitous Engineering (MUE '08), pp. 201–206, April 2008.
- [19] P. N. Belhumeur, J. P. Hespanha, and D. J. Kriegman, "Eigenfaces vs. fisherfaces: recognition using class specific linear projection," *IEEE Transactions on Pattern Analysis and Machine Intelligence*, vol. 19, no. 7, pp. 711–720, 1997.




















Submit your manuscripts at http://www.hindawi.com











Journal of Discrete Mathematics

