


**SUPPLEMENTAL FIGURE LEGEND (color reproduction only on the Web)**

**Figure S1.**



**Figure S1.** A) Genital system of Swiss male mice. B) CE duct segment of approximately 1.0cm length was isolated C) Segment was attached to FORT10 isometric force transducer (WPI, USA) connected to Transbridge 4M Transducer Amplifier (WPI, USA) to record the isometric tension development in a PC-based System (MP100, Biopac System INC., USA) and D) Data were analyzed off-line using the software AcqKnowledge version 3.5.7 (Biopac System INC., USA).